

Specifications of L-0187

●MICROSCOPE	
Type:	Galilean-Type Binocular Microscope
Magnification Changer:	3 Steps by drum rotation
Eyepieces:	15X
Total Magnifications:	10, 16, 25X
Field of View:	28.5, 14, 7mm
Pupillary Adjustment:	48-78mm
Focus Distance:	100 mm
Diopter Adjustment:	±6D
●SLIT ILLUMINATION	
Slit Width:	0-10mm
Aperture Diameter:	0.2, 1, 3, 4, 6, 10mm
Slit Length:	0-10mm
Slit Rotation:	±90°
Filters:	Cobalt blue, Heat absorption, Red-free, 1/2ND
Slit Inclination:	5° 10° 15° 20°
Lamp:	12V/50W Halogen Lamp(L0160H1)
●BASE	
Vertical Movement:	30mm
Longitudinal Movement:	70mm
Lateral Movement:	115mm
Fine Base Movement:	±10mm
●CHINREST	
Vertical Movement:	75mm
Fixation Target:	Tungsten lamp 6V 0.15A
Input Voltage:	AC. 100, 110, 220, 240V 50/60Hz
Power Consumption:	50VA
●DIMENSION & WEIGHT	
Weight:	20.5kgs
Dimension:	550(W)x350(D)x500(H)



Inami & Co., Ltd.

L-0187

3 STEPS MAGNIFICATION SLIT LAMP



L-0189
5 STEPS MAGNIFICATION SLIT LAMP



●Specifications subject to change without prior notice.

3 STEPS MAGNIFICATION SLIT LAMP L-0187



Inami & Co., Ltd.

BIJYOGI RESEARCH INSTITUTE
No. 3-14-13, Bijyogi, Toda-City,
Saitama-Pref, 335-0031, Japan
Tel. 81-48-422-0173 Fax. 81-48-421-3838
Website :http://www.inami.co.jp
E-Mail : ibri@inami.co.jp
®TRADE MARK OF INAMI & CO., LTD.

HEAD OFFICE
No. 24-2, Hongo 3-Chome,
Bunkyo-ku, Tokyo 113-0033, Japan
Tel. 81-3-3814-5916 Fax. 81-3-5684-2126
Website :http://www.inami.co.jp
E-Mail : corpinfo@inami.co.jp
®TRADE MARK OF INAMI & CO., LTD.

The L-0187 3 STEPS MAGNIFICATION SLIT LAMP provides practitioners with the utmost satisfaction thanks to its excellent image clearness and fine design.



MS-51
3 STEPS MAGNIFICATION SYSTEM
(10, 16, 25X, OBJECTIVE LENS f=100mm)



OCL-009
EYE PIECES 15X




YF-57
SPECIAL ADAPTOR BUILT-IN
YELLOW FILTER



VARIOUS KINDS OF BUILT-IN FILTERS:
COBALT BLUE, RED-FREE, 50%ND AND
HEAT ABSORPTION FILTERS



L-0542-C
ASSISTANT BINOCULAR MICROSCOPE



L-0541
BEAM SPLITTER



L-0544
TV CAMERA ADAPTOR(f=100mm)

