

## **Cannulae**

Re-usable

Disposable

AC

Lacrimal

Anaesthesia

Aspirating

Cystotomes

Capsule Polishers

Hydrodissectors

# Cannulae

## Retrobulbar Needles, 26g, box 10

RE-USABLE

0108008

50mm

0108007

38mm



## Rycroft Cannula, box 10

0108003

30g

0108004

27g



## Air Injection Cannula

0108010

27g, angled

0108091

27g, Frankfurt, curved

0108094

23g, curved

0108093

25g, standard, curved



## Anterior Chamber Irrigator

RE-USABLE

**0108088**

Southampton, 19g cannula, long

**0108084**

Bishop-Harmon Set, 19g cannula, connector and bulb

**0108083**

Bishop-Harmon, 19g cannula

**0108085**

Silicone Bulb

**0108087**

Connector



## Irrigating Vectis single hole, straight, 24g

**0108067**



## McIntyre-Binkhorst Cannula, blunt tip, 26g

**0108025**

right

\*(not shown)

**0108026**

left



# Cannulae

## Jensen Capsule Polisher

RE-USABLE

0108096

27g

0108097

23g



## Simcoe I/A Cannula, 23g

0108222

0.3mm aspiration port (shown)

0108223

reversed, 0.3mm aspiration port

0108224

0.4mm aspiration port

0108225


reversed, 0.4mm aspiration port

0108226

polishing tip, 0.3mm aspiration port

0108228

polishing tip, 0.4mm aspiration port



## Lacrimal Cannula, bulbous end, 19g

0108033

short

0108032

long



## Lacrimal Cannula, reinforced, 23g

RE-USABLE

0108035

curved

0108037

straight



---

## Lacrimal Cannula

0108040

Duggan, curved, 23g

0108036

Harrison Butler, 24g

0108038

straight, 23g



---

## Cannula Holder male/female

0108089




## Anaesthesia Needles

DISPOSABLE

**0180017**  
Peribulbar Atkinson  
26g x 19mm

BOX OF    DETAIL

10




**0181275**  
Retrobulbar Sharp  
25g x 38mm

10



**0181637**  
Retrobulbar Atkinson  
25g x 38mm

10




**0181638**  
Retrobulbar Atkinson  
23g x 38mm

10



**0185027**  
Peribulbar Atkinson  
25g x 22mm

10




**0185108**  
Peribulbar/Retrobulbar Atkinson  
25g x 32mm

10



**0185109**  
Peribulbar/Retrobulbar Atkinson  
27g x 32mm

10



**0185144**  
Peribulbar Atkinson  
23g x 22mm

10



**0185146**  
Facial Infiltration sharp  
30g x 16mm

10



**0185176**  
Sub Tenon curved  
19g x 25mm

10



DISPOSABLE

## Rycroft Anterior Chamber Cannula

**0181273**  
30g x 22mm  
angle 5mm from end

**0181274**  
27g x 22mm,  
angle 5mm from end

**0181639**  
25g x 22mm  
angle 7mm from end

**0185025**  
23g x 25mm  
angle 7mm from end

**0185175**  
30g x 22mm  
straight

**0185006**  
27g x 22mm  
straight

BOX OF DETAIL

10



10



10




10



10



10



## Aspirating Cannula

**0180111**  
Binkhorst, left hand  
25g x 22mm, 'U' shape

**0180113**  
Binkhorst, right hand  
25g x 22mm, 'U' shape

**0185090**  
Viscoexpression  
25g x 22mm, curved Corydon

**0185191**  
Hydroexpression  
27g x 22mm, curved Corydon

10




10



10



10



# Cannulae

## Capsule Polishers

DISPOSABLE

**0181604**

0.3mm top port, 21g x 22mm,  
sandblasted tip, 35 degree angle  
9mm from end

**0181605**

as above 23g

**0185005**

Silicone Tip, 27g x 22mm,  
45 degree angle 7mm from end

**0185041**

Silicone Tip, 23g x 22mm,  
45 degree angle 7mm from end

**0185121**

Silicone Tip, 23g x 22mm,  
45 degree angle 9mm from end

BOX OF DETAIL

10




10



10



10



## Hydrodissectors

**0181592**

Flattened, 30g x 22mm,  
angle 7mm from end

**0185029**

'J' Shaped,  
27g x 22mm

**0185037**

Flattened 25g x 22mm,  
angled 8mm from end


**0185155**

Flattened horizontally,  
27g x 22mm, Helsinki, curved

10






10



10



10





## Hydrodissectors (con't)

DISPOSABLE

**0185099**  
 Flattened vertically,  
 27g x 22mm, curved, Sauter


BOX OF    DETAIL

10




**0185107**  
 Hydrodissector/Hydrodelineator,  
 25g x 22mm, tapered,  
 angle 8mm from end

10



**0185158**  
 Cortical Cleaving,  
 27g x 22mm, flattened,  
 angle 11mm from end


10



## I/A Cannulae

**0185047**  
 Simcoe, 23g x 16mm,  
 0.3mm aspiration port

5



**0185048**  
 Simcoe, 23g x 12mm,  
 0.3mm aspiration port

5

\*(not shown)

**0185222**  
 Bi-manual, 23g x 17mm,  
 curved 45 degree,  
 0.4mm aspiration port

10



**0185223**  
 Bi-manual, 21g x 17mm,  
 curved 45 degree,  
 0.4mm irrigation port

10



**0181611**  
 21g x 38mm, formed

10

\*(not shown)

# Cannulae


## Irrigating Cystotomes

DISPOSABLE

**0181610**  
25g x 16mm, formed



BOX OF    DETAIL

10



**0181602**  
Side-cutting, 25g x 16mm, formed

10





**0181618**  
27g x 16mm, formed

10




**0181631**  
27g x 16mm, straight

10




**0185030**  
Intercapsular/Endocapsular, 25g x 16mm, curved

10



**0185042**  
27g x 16mm, formed

10



## Irrigating Vectis

**0181619**  
25g, single port

5



**0181626**  
25g, double ports

5



DISPOSABLE

## Lacrimal Cannula

**0181276**

26g x 28mm, curved

**0181615**

21g x 35mm, curved (not shown)

BOX OF DETAIL

10



10



## Anterior Chamber Maintainers

**0185061**

self-retaining,  
20g x 3.5mm


**0185163**

micro-etched,  
20g x 4mm

3



10



## Subretinal Fluid Cannulas

**0185021**

20g, straight, silicone tip, 6mm  
beyond end, total length 35mm

**0185035**

20g, curved, silicone tip, 6mm  
beyond end, total length 35mm

5



5



## Irrigating Dialer

**0181286**

30g x 22mm, angle 10mm  
from end

10



## Microsurgical Handle

**0181012**

luer male/female

20

