



HS JOHN WEISS
INTERNATIONAL

Disposables & Drapes

Eye Shields

Cauteries

Sponges

Drapes

Eye Shields, sterile, box 100

EYE SHIELDS

0111101

Universal, adult

0111104

right

0111105

left



*(actual size not shown)

Eye Patches, black silk with cord (10), non sterile

0120300



Cautery

Weiss Disposable, sterile, box 10

0111122



CAUTERIES & EPILATION NEEDLES


Cautery Burner Weiss

0120104



Cautery Burner Dixey

0120105



Cautery Burner Red Heat

0120108



Epilation Needles, straight, box of 50, sterile

0120133



Sugi®Swabs

SPONGE

0111001
sterile, qty 20 x 10


0111002
non sterile, qty 500



Sugi®Spears and Strips

0111004
sterile, qty 20 x 10

0111008
sterile, qty 20 x 10



Weiss Cellulose Spears

1701002
sterile, qty 20 x 5


1701003
sterile, qty 20 x 10



Weiss PVA Spears

1701012
sterile, qty 20 x 5

1701013
sterile, qty 20 x 10




Weiss Instrument Wipe

SPONGE

1701010


sterile, qty 20,
75mm x 75mm



Weiss Diamond Knife Cleaning Blocks

1701009

sterile, qty 10



Weiss Eye Wicks

1701006

sterile, qty 10 X 2


Weiss Eye Wick Drains

1701004

80cc bag, sterile, qty 10

1701005

400cc bag, sterile, qty 10



John Weiss Ophthalmic Drapes

- Super soft non-woven fabric - drapes like real cotton
- Breathable, lightweight and moisture repellent for maximum patient comfort
- Adhesive material is odourless and designed for atraumatic removal
- Generous dimensions
- Easy-to-follow instructions
- Channel design bag for accurate fluid collection, made from non-reflective translucent polymer

Ophthalmic Drape Full Body, Incise Non-Woven, Box of 10

1601087


Overall size: W160 x L250cm

Incise area: W10 x L8cm

Dual channel bag: W53 x L30cm

Tenting Band: 1

Wick: 1



Ophthalmic Drape Large, Incise Non-Woven, Box of 10


1601084

Overall size: W160 x L160cm

Incise area: W10 x L8cm

Channel bag: W34.5 x L25.5cm

Wick: 1



DRAPES

Ophthalmic Drape Medium, Incise Non-Woven, Box of 10


1601072

Overall size: W130 x L115cm

Incise area: W10 x L8cm

Fluid bag: W20 x L30cm

Wick: 1



Ophthalmic Drape, Medium Plastic, Box of 10


1601071

Overall size: W100 x L127cm

Incise area: W11.5 x L9cm

Fluid bag: W53 x L30cm

Wick: 1



Ophthalmic Drape, Mini Plastic, Box of 10

1601095

Overall size: W40 x L40cm

Incise area: W11.5 x L9cm

