

Oculoplastic & Lacrimal

Probes & Dilators
DCR & Intubation Sets
Jones Tear Duct Tubes
DCR Instruments
Orbital Instruments
Enucleation &
Evisceration
Spheres & Shells

Liebreichs (Bowmans) Lacrimal Probes, silver

PROBES & DILATORS

0105008

set of four

0105016

3/4

0105013

1/2

0105012

00/0

0105011

0000/000



Pigtail Probe with holes, hexagonal handle

0105010



Nettleship Lacrimal Dilator

0105033

single ended



0105034

double ended



Wilder Lacrimal Dilator

0105036

fine taper

20mm length



0105037

medium taper

16mm length



Lacrimal Dilator, infant

PROBES & DILATORS

0105038

delicate taper
22mm length



Castroviejo Lacrimal Dilator, double ended

0105040

22mm/18mm
taper length



Hays Punctum Finder

0105154




DCR Bodkins Set, angled, box 2, sterile

DCR & INTUBATION SETS


0111026

angled



0111028


straight



Weiss Lacrimal Intubation Set, 18cm probes, box 2, sterile


DCR & INTUBATION SETS

0111027



Crawford Lacrimal Intubation Set, olive shaped probe tips, box 2, sterile

0111029



0111030

Crawford Retrieval Device





JONES TEAR DUCT TUBES

2.2mm Pyrex Tubing		Flange Diameter & Tube Length mm
0110019	Bevelled	3.0 x 10
0110029	Bevelled	3.0 x 11
0110020	Bevelled	3.0 x 12
0110030	Bevelled	3.0 x 13
0110009	Bevelled	3.0 x 14
0110031	Bevelled	3.0 x 15
0110022	Bevelled	3.0 x 16
0110003	Bevelled	3.5 x 10
0110014	Bevelled	3.5 x 11
0110004	Bevelled	3.5 x 12
0110005	Bevelled	3.5 x 13
0110006	Bevelled	3.5 x 14
0110007	Bevelled	3.5 x 15
0110008	Bevelled	3.5 x 16
0110018	Bevelled	3.5 x 17
0110026	Bevelled	3.5 x 18
0110027	Bevelled	3.5 x 19
0110032	Bevelled	4.0 x 12
0110033	Bevelled	4.0 x 13
0110034	Bevelled	4.0 x 14
0110010	Bevelled	4.0 x 15
0110011	Bevelled	4.0 x 16
0110039	Bevelled	4.0 x 17
0110012	Bevelled	4.0 x 18
0110040	Bevelled	4.0 x 19
0110041	Bevelled	4.0 x 20
0110054	Bevelled	4.0 x 21
0110043	Bevelled	4.0 x 22
0110059	Angled & Bevelled	3.0 x 12
0110060	Angled & Bevelled	3.0 x 16
0110048	Angled & Bevelled	3.0 x 17
0110058	Angled & Bevelled	3.5 x 12
0110016	Angled & Bevelled	3.5 x 14
0110017	Angled & Bevelled	3.5 x 15
0110055	Angled & Bevelled	3.5 x 16
0110049	Angled & Bevelled	3.5 x 17
0110062	Angled & Bevelled	3.5 x 19
0110024	Angled & Bevelled	4.0 x 14
0110025	Angled & Bevelled	4.0 x 15
0110050	Angled & Bevelled	4.0 x 16
0110051	Angled & Bevelled	4.0 x 17
0110053	Angled & Bevelled	4.0 x 18
0110052	Angled & Bevelled	4.0 x 19
0110042	Straight	4.0 x 21
0110044	Straight	4.0 x 24

* Jones Tear Duct Tubes also available in frosted version. Please specify when ordering.

* Tubes can be made to specification and with suture holes. Please specify flange diameter and length.

JONES TEAR DUCT TUBES



Other Products		Flange Diameter & Tube Length mm	
0110036	Bevelled Poly Tube	4.5 x 18	
0110037	Tubing	2.2 x 3.5 x 30mm	
0110045	Nylon Cleaning Rod		
0110046	Jones Gold Dilator No. 1		
0110047	Jones Gold Dilator No. 2		
0112052	Storage Tray		
0112055	Sterilising Tray (shown)		
0112056	Measuring Tray (shown)		

Sterilising Tray



Measuring Scale

SETS

0110002 STANDARD SET

x 2	3.0 x 14		
x 2	3.5 x 12		
x 2	3.5 x 13		
x 2	3.5 x 14		
x 2	3.5 x 15		
x 2	4.0 x 15		
x 1	4.0 x 16		
x 1	4.0 x 18		
x 2	Jones Poly Tubes (4.5 flange)		
x 5	Nylon Cleaning Rods		
x 1	Jones Gold Dilator 1		
x 1	Jones Gold Dilator 2		

0110013 ASSORTED SET

x 1	3.0 x 10		
x 1	3.0 x 12		
x 1	3.0 x 14		
x 1	3.0 x 16		
x 1	3.5 x 11		
x 1	3.5 x 13		
x 1	3.5 x 15		
x 1	3.5 x 17		
x 1	4.0 x 15		
x 1	4.0 x 14 angled		
x 1	4.0 x 15 angled		
x 1	3.5 x 14 angled		
x 1	3.5 x 15 angled		
x 2	Jones Poly Tubes (4.5 flange)		
x 5	Nylon Cleaning Rods		
x 1	Jones Gold Dilator 1		
x 1	Jones Gold Dilator 2		

SETS (con't)

JONES TEAR DUCT TUBES

0110038 SPECIAL SET	
x 1	3.0 x 10
x 1	3.0 x 11
x 1	3.0 x 12
x 1	3.0 x 13
x 1	3.0 x 14
x 1	3.0 x 15
x 1	3.5 x 9
x 1	3.5 x 11
x 1	3.5 x 17
x 1	3.5 x 18
x 1	4.0 x 12
x 1	4.0 x 13
x 1	4.0 x 14
x 1	4.0 x 17
x 1	4.0 x 19
x 2	Jones Poly Tubes (4.5 flange)
x 5	Nylon Cleaning Rods
x 1	Jones Gold Dilator No. 1
x 1	Jones Gold Dilator No. 2

Weiss Sphenoidal Punch, fine

DCR INSTRUMENTS

0101487
upward cut, 2mm

0101488
downward cut, 2mm

0101489
upward cut, 3mm

0101490
downward cut, 3mm

0101495
curved



*(actual size not shown)

Kerrison Punch

DCR INSTRUMENTS

0204149

4mm

*not shown

Killian Nasal Speculum

0206020



Friedman Bone Nibbler Forceps, curved

0120308



Zoellner Suction Tube, for DCR

DCR INSTRUMENTS

0120311



Knapp Lacrimal Sac Retractor, four blunt prongs

0107052



Rollet Lacrimal Sac Retractor, four sharp prongs

0107059



Stevenson Lacrimal Sac Retractor

0107057



Rollet Rugine

DCR INSTRUMENTS

0114029



Rollet Rugine and Periosteal Elevator, double ended instrument

0114299



Freer Elevator, double ended, sharp/blunt, 182mm length

0114262



*(actual size not shown)

Traquair Periosteal Elevator, for DCR

0114263



Howarth Periosteal Elevator, for DCR, 215mm length

0114264



*(actual size not shown)

Stallard Dissector, blunt for DCR

0114162



Mallet 210mm length

DCR INSTRUMENTS

0114312



*(actual size not shown)

Straight Chisels

0114084

3mm

0114085

4mm

0114086

5mm

0114087

6mm

0114088

7mm



*(actual length 180mm)

Curved Chisels

0114103

3mm

0114104

4mm

0114105

5mm

0114106

6mm

0114107

7mm



*(actual length 180mm)

Osteotomes

DCR INSTRUMENTS

0114143

3mm

0114144

4mm

0114145

5mm

0114146

6mm

0114147

7mm

*(actual length 180mm)



Watson Orbital Retractors, Maleable, set of four

ORBITAL INSTRUMENTS

0107072



*(actual size not shown)

Fison Retractor, for detachment

0107092



Downs Rootman Orbital Surgery Set

ORBITAL INSTRUMENTS

0107112

SET QUANTITIES

- 0107132
Forceps, Tumour Holding, square notched jaws
- 0107131
Forceps, Tumour Holding, round cup jaws
- 0107130
Hook, Tumour Holding, blunt
- 0107129
Hook, Tumour Holding, sharp
- 0107128
Retractor, Vessel
- 0107127
Dissector/Elevator, curved 5mm
- 0107126
Dissector/Elevator, curved 4mm
- 0107125
Dissector/Elevator, curved 3mm
- 0107124
Dissector/Elevator, straight 4mm
- 0107123
Retractor, Temporal Fossa, right
- 0107122
Retractor, Temporal Fossa, left
- 0107121
Elevator, Periosteal, round 5mm
- 0107120
Elevator, Periosteal, round 7.5mm
- 0107119
Elevator, Periosteal, double ended
- 0107118
Retractor, skin
- 0107117
Retractor, 15mm wide curved general
- 0107116
Retractor, 8mm wide curved general
- 0107115
Retractor, 32mm deep flared end
- 0107114
Retractor, 16mm deep
- 0107113
Retractor, 10mm deep

1

1

1

1

1

1

1

1

1

1

1

1

1

4

2

2

1

2

2



ALSO AVAILABLE INDIVIDUALLY

(suggested tray 0112004)

*(actual size not shown)

Trelat Enucleation Spatula

ENUCLEATION & EVISCERATION

0106045



Mule Evisceration Scoop

0114240



Elliott Evisceration Scoop

0114241



Moorfields Globe Introducer

0120303



Foster Enucleation Snare

0114325

0114326

snare wire 29 gauge
(2 metres)

0114327

snare wire 26 gauge
(2 metres)



*(actual size not shown)

Castroviejo Implant Acrylic

SPHERES & SHELLS

0120072
regular
11.5mm depth

0120073
deep
13mm depth



Implant Retaining Shell

0120074



Allan Implant Acrylic, non magnetic

0120077
19mm

0120078
21mm



Mule Globe Vitreous

0120042
12mm

0120043
14mm

0120044
16mm

0120045
18mm

0120046
20mm



Wright Globe Acrylic

SPHERES & SHELLS

0120047

12mm

0120048

14mm

0120049

16mm

0120050

18mm

0120051


20mm

0120053

21mm

0120052

22mm



Conformers perforated, plain glass, pack of 5

0120013

right, small

0120020

right, medium

0120022

right, large

0120019

left, small

0120021

left, medium

0120023

left, large



Conformers plain glass, pack of 5

0120018

right, small

0120015

right, medium

0120017

right, large

0120012

left, small

0120014

left, medium

0120016

left, large

