


HS JOHN WEISS
INTERNATIONAL



Vitreoretinal

Forceps

Scissors

Membrane
Instruments

Subretinal and
ILM Instruments

Vitreoretinal

Ingenuity You Inspire

The range of vitreoretinal micro forceps and scissors from John Weiss offers exceptional quality in a wide variety of instrument styles and sizes. Each instrument is handcrafted by our expert instrument makers from the highest grade materials, resulting in consistent accuracy, outstanding precision, exquisite feel and movement.

The Weiss style rounded grasp area of the delicate titanium handle enables easy 360° rotation and extremely smooth movement for high precision surgical manoeuvres. A wide variety of instrument tips are available to meet every surgical need. Tips are made of a special stainless steel, which is highly corrosion resistant, resulting in smoother operation and fewer repairs.

The Weiss collection is available in 20 gauge and in 23 and 25 gauge trocar compatible versions.



MICRO SCISSORS

Straight Scissors

Sharp tipped scissors with long blades

2003002


20g 3mm

2003038

23g trocar 3mm

2003044

25g trocar 1.25mm



Curved Scissors

Sharp tipped scissors with long blades

2003003


20g 3mm

2003039

23g trocar 1.25mm

2003045

25g trocar 1.25mm



135° Angled Scissors

Sharp tipped horizontal style scissors angled 135°, 2mm from tip

2003004

20g



125° Angled Scissors

Sharp tipped horizontal style scissors angled 125°, 2mm from tip

2003006

20g



60° Angled Scissors

Sharp tipped fine vertical style scissors with long blades

2003007

20g 1.9mm



MICRO SCISSORS

45° Angled Scissors

Sharp tipped fine vertical style scissors with long blades

2003005


20g 1.75mm

2003040

23g trocar 0.6mm

2003046

25g trocar 0.5mm




*2003040 demonstrates how the width of a Trocar tip does not exceed the gauge of the instrument.

Pointed End Gripping Forceps

Fine pointed forceps with platform

2003061

23g trocar




MICRO FORCEPS

End Gripping Forceps

Blunt 4mm long jaws with small gripping platform at each tip

2003012

20g 4mm



MICRO FORCEPS

Micro End Gripping Forceps

Very fine blunt long jaws with small platform at each tip

2003013


20g 3.5mm

2003041

23g trocar 3mm

2003062

25g trocar 1.25mm




Angled Forceps

Blunt tipped forceps with small platforms, angled 135°, 2mm from the tip

2003015

20g




Pick Forceps

Straight long blunt jaws with large platform on each tip

2003016

20g 4mm



Serrated Forceps

Blunt tipped forceps with long jaws

2003017

20g 4mm

2003042

23g trocar 3mm

2003047

25g trocar 1.5mm




Foreign Body Forceps

Reverse action basket style forceps with 3 gripping arms

2003019

20g



MICRO FORCEPS

Long Tipped Pick Forceps

Small pick forceps with platform

2003018

20g 0.8mm



Offset Gripping Forceps

Platform is pointed to one side for grasping along the retinal surface

2003020

20g

2003063

23g trocar 3mm

2003049

25g trocar



ILM Pick

Micro pick for lifting or removal of ILM or other membranes

2003058

20g

2003068

23g trocar

2003057

25g trocar



MEMBRANE INSTRUMENTS

SUBRETINAL & ILM

Pick Forceps


Jaws with small pick and platform at the tip for lifting membranes

2003064

23g trocar 3mm

2003048

25g trocar 1.5mm



Long Angled Scissors

Sharp tipped horizontal style scissors with thin blades angled 135°, 3.5mm from tip

2003022

20g



Angled Long Gripping Forceps

Delicate forceps with jaws angled 135°, 3mm from tip

2003023

20g



Fine Gripping Forceps

Fine long jaws with small platform at each tip

2003024


20g 3.5mm

2003043

23g trocar

2003050

25g trocar



Ducournau Fine Gripping Forceps

Delicate forceps angled 135°, 3.5mm from tip

2003025

20g



Serrated ILM Forceps

Shortened serrated forceps with pick at the tip

2003065

23g trocar



SUBRETINAL & ILM

Naito ILM Forceps

Fine long jaws with a 0.7mm pick at the tip

2003026


20g 3.3mm

2003066

23g trocar 2.2mm

2003067

25g trocar 2.2mm




Delicate Grasping Forceps

Delicate long, straight forceps with fine platform at each tip

2003028

20g 3.5mm



Vitreoretinal Flushing System

2003069

