

Extra long jaws help to advance the IOL deep into the cartridge

Angled and straight IOL loading forceps

[www.johnweiss.com](http://www.johnweiss.com)

 **JOHN WEISS  
INTERNATIONAL**

### John Weiss IOL loading forceps offer some key benefits:

- Extra long jaws enable you to advance the IOL deep into the cartridge
- Smooth atraumatic rounded tips reduce risk of damage to the IOL
- Tips form a 'v' shape when closed, facilitating better engagement with the IOL when advancing
- Angled and straight options available in either stainless steel or titanium.

If you'd like to arrange to see them please contact us on **+44(0)1908 340000** or email **[sales@johnweiss.com](mailto:sales@johnweiss.com)**



**0109448** Titanium  
**0101610** Stainless Steel

IOL loading forceps, angled, 15mm tip to angle jaws, overall length 105mm



**0109449** Titanium  
**0101609** Stainless Steel

IOL loading forceps, straight, 15mm jaws, overall length 110mm



Please quote the product code when ordering